

Adaptaid⁺

HEPA

Ventilated

Cabinet for Sterile Containers



Cabinets designed for
the safe storage of
medical instruments.

Endoscopy

Sterilization - MDR

ENT

adaptaid.com



HEPA Ventilated Cabinet for Sterile Containers MEDICASE

Advantages

- + **Cabinet designed to add traceability functions** and ensure constant monitoring of equipment.
- + **Large door windows** for better visibility.
- + **Evenly distributed LED lighting.**
- + **Flexible storage system** for containers or bins.
- + **Quick access to filter change.**
- + **Locking device** to prevent intrusion and theft.
- + **Manual door opening device** in case of emergency.
- + **Height-adjustable shelves** to maximize space for different sizes of containers or bins.



Single or Double Door Cabinets

Models	Doors	Shelves	Dimensions inches W x D x H	Dimensions millimeters W x D x H
EDS-024-3027	1	6	30 x 27 x 72	762 x 686 x 1829
EDS-024-4227	2	6	42 x 27 x 72	1067 x 686 x 1829

Hepa Ventilated Sealed Cabinet with Controlled Openings to Ensure an Aseptic Environment for Storing Sterile Containers

Structure

- + 100% stainless steel

Shelves

- + Height adjustable shelves with ventilation openings
- + Pull-out for easy access (option)

Filters

- + HEPA filters 99.97%
- + Replacement indicator and pre-filter

Ventilation

- + Continuous positive cabinet pressure

Doors

- + Tempered glass
- + Completely sealed
- + Open door warning alarm
- + Key lock (Electronic lock on order)

Other Features

- + LED light

Compliance

- + Designed according to the standards in force and industry recommendations for medical environments.
- + Made in Canada



The addition of an ADAPTAID+ traceability system ensures constant monitoring of the cabinet

The ADAPTAID+ Traceability System ensures constant monitoring of the status of your endoscopes, while having control over their accessibility at all times. The Department Manager can access the logbook locally or remotely to see all recorded events. The system is equipped with high-performance tools such as a color touch screen, a barcode reader and an RFID reader.

Traceability System (on order)

- + Integrated system
- + 10 inch (254 mm) color touch screen
- + Bar code RFID reader
- + Events and alarms recorded in the logbook
- + Color codes and pictograms for easy identification
- + Recording of stored instruments
- + Easy Identification of instrument status
- + Alerts for filter changes
- + Controlled users access
- + Controlled storage time
- + Monitoring of cabinet maintenance
- + Electronic lock
- + Open door alert



Adaptaid⁺

📍 448 rue Bourque
Repentigny, Québec
J5Z 5A2 Canada

☎ 450-585-7272

@ info@adaptaid.com



THIS DOCUMENT IS THE PROPERTY OF ADAPTAIDE INC. ANY REPRODUCTION, USE OR COPY FOR THE PURPOSE OF MANUFACTURING OR SELLING ANY PRODUCT IS STRICTLY PROHIBITED WITHOUT WRITTEN CONSENT.

2023 ADAPTAIDE INC. ALL RIGHTS RESERVED