

Adaptaid⁺

Storage Cabinet for Endocavity Probes MEDIPROB



Ventilated cabinet designed to store transesophageal echocardiogram (TOE/TEE) probes efficiently and safely.

**Endoscopy
Operating Room**

adaptaid.com



MEDIPROB

Endocavity Probe Ventilated Storage Cabinet



Storage with Fixed Probe Holders

Models	Probes	Doors	Traceability	Dimensions W x D x H	
				inches	mm
EDS-028-1810EN	4	1		18 x 10 x 80	457 x 254 x 2032
EDS-028-2410EN	6	1		24 x 10 x 80	610 x 254 x 2032
EDS-028-3010EN	8	1		30 x 10 x 80	762 x 254 x 2032
EDS-028-3810EN	10	2		38 x 10 x 80	965 x 254 x 2032
EDS-028-4610EN	12	2		46 x 10 x 80	1168 x 254 x 2032
EDS-028-LT-1810EN	4	1	✓	18 x 10 x 80	457 x 254 x 2032
EDS-028-LT-2410EN	6	1	✓	24 x 10 x 80	610 x 254 x 2032
EDS-028-LT-3010EN	8	1	✓	30 x 10 x 80	762 x 254 x 2032
EDS-028-LT-3810EN	10	2	✓	38 x 10 x 80	965 x 254 x 2032
EDS-028-LT-4610EN	12	2	✓	46 x 10 x 80	1168 x 254 x 2032

Efficient and Secure Probe Storage

- + **Direct front access** to probes.
- + **Individual storage system.**
- + **Evenly distributed LED lighting.**
- + **Soft molded surfaces** to prevent damage to probes.
- + **Sealed HEPA filtered cabinet** to maintain positive pressure and ensure an aseptic environment.



Specifications

Frame

- + 100% stainless steel construction

Filters

- + HEPA filter 99.97% for the inner chamber
- + Air inlet and outlet pre-filters
- + Filter replacement indicator

Air Filtration

- + Continuous airflow using positive pressure
- + Vertical ventilation from top to bottom, evenly distributed
- + Ventilation blockage / failure indicator

All You Need to Ensure Safe Use of Endocavity Probes

Probe Holders

- + Four to twelve fixed holders for endocavity probes
- + Front access for filter change
- + Probe holders with soft probe rest

Doors

- + Available with left-hand or right-hand opening
- + Tempered glass
- + Perfectly sealed
- + Open door warning alarm
- + Mechanical lock

Additional Options

- + Traceability system
- + Code or RFID lock
- + Custom holders for special applications

Other Features


- + Wall anchor bar
- + LED lighting

Compliance

- + Designed according to the standards and recommendations in force in the medical industry.
- + Made in Canada



Adaptaid⁺

 448 rue Bourque
Repentigny, Québec
J5Z 5A2 Canada

 450-585-7272

 info@adaptaid.com



THIS DOCUMENT IS THE PROPERTY OF ADAPTAIDE INC. ANY REPRODUCTION, USE OR COPY FOR THE PURPOSE OF MANUFACTURING OR SELLING ANY PRODUCT IS STRICTLY PROHIBITED WITHOUT WRITTEN CONSENT.

2023 ADAPTAIDE INC. ALL RIGHTS RESERVED