

Adaptaid<sup>+</sup>

---

# ERGOFLEX

## Sterile Wrap Rack



Sterilization/MDRD

## Retractable Sterile Wrap Rack

Installed directly on the base of the ERGOFLEX preparation and packing table, this retractable sterile wrap rack optimizes the storage of sterile wraps by reducing the space required on work surfaces, walls and floors.



© ADAPTAID  
PATENT PENDING


### Solid Structure

Robust, torsion-resistant 304 stainless steel structure. Four cylindrical bars hold the sterile wraps in place with retaining rings.



#### OPEN STERILE WRAP RACK

- + Quick and effortless opening
- + Direct access to wraps
- + Depth of 19 7/8 inches (505 mm) when open



#### CLOSED STERILE WRAP RACK

- + Closes with a light push
- + Increased lateral clearance
- + Depth of 12 inches (305 mm) when closed



## An Ergonomic, User-Friendly and Productive Solution

The ERGOFLEX retractable sterile wrap rack provides fast, close and ergonomic access to sterile wraps. What's more, its proximity to the table increases the comfort and efficiency of the sterile processing staff in carrying out the packaging process of medical devices. Among its many benefits are:

- + Ergonomic organization of sterile wraps.
- + Increase productivity by enabling a continuous workflow.
- + Optimization of available floor space.
- + Reduced musculoskeletal injuries thanks to lateral access to sterile wraps.
- + Reduced footprint thanks to the retractable frame.
- + Compatible with various sizes of sterile wraps (up to 48" wide (1219 mm)).

### Additional Options

- + Telescopic extension bar 40.4 to 51.5 inches (1026 to 1308 mm)
- + Fixed, non-retractable frame (model RDM-010-SCF)



MODELS	DIMENSIONS			*SPACE SAVING	
	INCHES	W X D X H	MM	INCHES	MM
RDM-010-SC	39.6 x 19.8 x 45.4	1006 x 505 x 1228		7.8	198
RDM-010-SCF	39.6 x 20 x 45.4	1006 x 508 x 1228		s/o	s/o

\* Lateral space gained when the retractable rack is closed.

# Adaptaid<sup>+</sup>

High-Performance Medical Equipment  
Designed for Healthcare Professionals

**For more information about our products, contact our customer service department.**



CONTACT US  
+1 450-585-7272  
customerservice@adaptaid.com



©2026 ADAPTAIDE INC. ALL RIGHTS RESERVED

[adaptaid.com](https://adaptaid.com)