

Adaptaid⁺

Endoscope and Probe Transport Tray



Sterilization/MDR
Ultrasound

Endoscopy
Operating Room

Sturdy and Durable Transport Tray

Designed specifically for the medical industry, this container with a transparent lid is equipped with two side handles for easy handling and a security seal device. A reusable label (clean/soiled) can be inserted under the lid to indicate the status of the instrument inside.



Compatible with Washer-Disinfectors

All components of this tray can be cleaned in washer-disinfectors at temperatures up to 115°C.

Tray Features

Compatible with large flexible probes and endoscopes.

| | |
|--------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Structure | Translucent gray polycarbonate molded tray with two side handles and a security seal device. Robust, durable, and eco-friendly, it requires no consumables. |
| Lid | Molded from clear polycarbonate designed to provide a quick view of the tray's contents. |
| Dimensions | W x D x H: 22.9 x 18 x 5 inches, 582 x 457 x 127 mm |
| Cleaning | Easy to clean and compatible with washer-disinfectors. |
| Accessories | Reversible and reusable identification label (clean/soiled). |



Eco-Responsible

The concept of the tray was developed to eliminate the use of consumables, so each component is cleanable and reusable. The use of this resealable tray with a reusable lid requires no disposable protective film. This advantage reduces waste and promotes sustainable and more economical practices.



Solutions Compatible with Transport Trays

- + Transport carts for 5 or 6 trays
- + MEDISTOR endoscope cabinets
- + Pass-through cabinets for trays
- + MEDIPROB cabinets



Designed for Transporting Flexible Endoscopes and Probes

This container facilitates the transport of probes as well as flexible endoscopes, such as gastroscopes, colonoscopes, bronchoscopes, cystoscopes, and enteroscopes.



Adaptaid⁺

High-Performance Medical Equipment
Designed for Healthcare Professionals

For more information about our products, contact our customer service department.



CONTACT US
+1 450-585-7272
customerservice@adaptaid.com



©2026 ADAPTAIDE INC. ALL RIGHTS RESERVED

adaptaid.com