

Adaptaid⁺

Dialysis Wall Box Station



Dialysis
Medical Water

Dialysis Wall Box Station

This stainless steel wall-mounted station has been redesigned down to the smallest detail to provide a reliable and durable solution capable of performing one or two treatments, depending on the configuration selected.

All of the station's water distribution components are designed to maintain a clean and safe environment, offering hospitals a premier solution for obtaining the medical-grade water required for hemodialysis treatments.



High-Performance and Safe

Thanks to the separation of purified water supply and wastewater drainage connections, this station provides filtered water that helps protect patients from infections and staff from potential contamination.



A Long-Lasting Solution with Many Benefits

- + A station with single or dual supply connections that allows for the treatment of one or two patients at a time
- + A versatile solution compatible with a wide variety of water and connection options, including high-purity water connections and raw water connections (optional).
- + Replacement-compatible housing: its dimensions fit easily into the openings of the recessed wall boxes being replaced.
- + Minimizes the risk of cross-contamination thanks to the wastewater drainage compartment, which prevents splashing onto the pure water connections and other components.
- + Corrosion-resistant materials for a long service life and rock-solid reliability.
- + Removable and easily cleanable wastewater drainage compartment.
- + Smooth, easy-to-clean surfaces.

Features of Wall Box Stations

High-Quality Components That Ensure Performance and Durability

Materials	16-gauge 316 stainless steel housing and 316 stainless steel fittings
Pure water inlet	Choice of single or double inlet tees for 25 mm standard PEX tubing (rear)
Pure water outlet	Single or double Swagelok QC4 female quick-connect fittings, depending on the model selected
Drain connection	Single or double Swagelok QC6 male quick-connect fittings, depending on the model selected
Dimensions	Internal dimensions (W x H): 12.5 x 13.5 inches (318 x 343 mm) External dimensions (W x H): 14.875 x 15.75 inches (378 x 400 mm) Wall-mounting frame: 1.375 inches (35 mm) wide
Installation	The installation kit includes mounting hardware for standard wall boxes.

Available Configurations

- WTR-042-P001 – Single 25 mm inlet tee, single outlet drain port
- WTR-042-P002 – Double 25 mm inlet tee, double outlet drain port
- WTR-042-P003 – Single 20 mm inlet tee, single outlet drain port
- WTR-042-P004 – Double 20 mm inlet tee, double outlet drain port
- WTR-042-P005 – Single inlet tee 32 mm, single outlet drain port
- WTR-042-P006 – Double inlet tee 32 mm, double outlet drain port
- WTR-042-P009 – Single blank tee with drain block
- WTR-042-P010 – Double blank tee with drain block
- WTR-042-P011 – Single raw water port, dual outlet drain port
- WTR-042-P012 – Double raw water port, dual outlet drain port

* The availability of certain configurations may vary.



Adaptaid⁺

High-Performance Medical Equipment
Designed for Healthcare Professionals

For more information about our products, contact our customer service department.



CONTACT US
+1 450-585-7272
customerservice@adaptaid.com



Proudly made in
CANADA

©2026 ADAPTAIDE INC. ALL RIGHTS RESERVED

adaptaid.com